[www.arc4health.com](http://www.arc4health.com)

**Applied Microcurrent Technology Limited**  
Unit 7, Woodland Enterprise Centre, Hastings Road  
Flimwell, East Sussex TN5 7PR, UK

© 2018 Applied Microcurrent Technology Limited  
Patent Numbers GB 2502133 & GB 2534681  
Ver. 4

## User manual

The complete pain management  
and tissue repair system

## What is the Arc4Health Device?

The Device is a wearable, non-invasive and drug-free treatment device that uses microcurrent therapy to support the management of pain and aims to repair tissue. The Device has EU Certification as a Class IIa Medical Device which means it can be used on humans. The Device is supplied subject to our terms and conditions of purchase available from our website. Please ensure that you read these prior to using the Device.

## Contraindications

- Do not use the Device if you have a latex or rubber allergy as some of the component parts and consumables may contain latex or rubber.
- Do not use the Device on people with cardiac pacemakers, implanted defibrillators, or any other implanted electronic devices, as this may cause electrical interference.
- Do not use the Device on babies, infants or children.
- The safety of microcurrent therapy during pregnancy has not been established, therefore we do not recommend it.
- The safety of microcurrent therapy with cancer has not been established. We do not recommend it for people with suspected or diagnosed cancer of any type.
- Users with suspected or diagnosed epilepsy should follow precautions recommended by physicians.
- Users with suspected or diagnosed heart disease and heart conditions of any type should follow precautions recommended by their physicians.
- Do not use the Device at the same time as a TENS machine.

## What's in the box?

- The Device
- Charger
- Arm Strap
- Conductive Pads (x4)\*
- User Manual

## Programmes


We always recommend you first run through the following 6-week treatment protocol in order to maximise ATP (adenosine triphosphate) production. Once the initial 6-week protocol finishes you can select the individual programme that best suits your condition (please see below for a description of each programme). Every few months, we recommend running through the 6-week protocol again in order to maintain optimum results (please note that this is a manufacturer's guidance and not a clinical recommendation).

As a general rule<sup>1</sup> : (i) acute phases of pain are defined as lasting less than 12 weeks; and (ii) chronic phases of pain are defined as lasting more than 12 weeks.

 Indicator lights

 On/Off and Programme button

 Play button

## 6 Week Protocol:

**P1** = Pain Management (Acute Phase) daily use for 2 weeks

**P2** = Pain Management (Chronic Phase) daily use for 1 week

**P3** = Tissue Repair (Acute Phase) daily use for 2 weeks

**P4** = Tissue Repair (Chronic Phase) daily use for 1 week

**Each programme lasts for 3 hours.**

<sup>1</sup>[www.britishpainsociety.org/media-resources](http://www.britishpainsociety.org/media-resources)

## Programming the Device

- On/Off button. Hold for 3 seconds to turn on or to turn off.
- All four indicator lights will very quickly show green, then red to indicate that your Device is functioning correctly.
- After this, lights will flash green corresponding with battery charge: 1 light = 25%, 2 = 50%, 3 = 75%, 4 = 100%
- After 3 seconds all lights will go out. Immediately after this, an amber light will indicate which programme was last used. The light will be green if you already have a connection. 1st light = **P1**, 2nd = **P2**, 3rd = **P3**, 4th = **P4**
- Press the left (On/Off) button, quickly, to skip through and select the desired programme from **P1** to **P4**.
- Press the right play button, quickly, to activate the selected programme. The unit will beep on activation of the selected programme.
- Once activated, the right play button will turn either green (connected) or amber (not connected). After 30 seconds the light will go out to conserve battery, but the programme indicator will continue to flash throughout the delivery cycle.

## How to put the Device on

**Only apply the Device to normal, intact, unbroken, clean, healthy skin.**

**Conductive Pads:** Fix to the arm strap using the press studs. Please ensure that the conductive pads are placed back onto the Perspex film after use to help prolong the adhesiveness.

**Arm Strap:** Attach the arm strap with the Device to the outside of the upper arm. Rotate the strap between arms each day. Please do not over tighten - the conductive pads are very sensitive only light pressure is required to maintain connectivity. After putting the Device into the arm strap's pouch, please ensure that the electrical lead in the arm strap is inserted fully into the bottom of the Device.

**Never strap the Device to an injured arm and where possible, alternate the arm that the Device is attached to on a daily basis and reposition the Device too\*\*.**

**\*Please Note:** Our new kits contain conductive pads which do not require conductive gel. Older kits may contain carbon rubber pads (x2) and a bottle of conductive gel (250ml). The carbon rubber pads will require the application of a small amount of conductive gel to ensure a good electrical connection with the skin. The conductive gel is propylene glycol free. Once the carbon rubber pads and conductive gel are phased out from the new kits and are replaced by the conductive pads, they will still be able to be purchased separately from our website.

Do not use any other gel in conjunction with the Device. Even some medical grade conductive gels can contain propylene glycol, which may react with the skin.

## Charging

Only use the charger supplied, or an original replacement charger ordered directly from Device Manufacturer.

Please charge the Device fully when you first receive the Device. One 4-5 hour charge can give up to 30 hours of use, enough for 10 delivery cycles.

Insert the charger lead into the base of the Device. Plug into power socket. Two beeps indicate charging.

Battery charge is indicated by the corresponding number of flashing green 'indicator lights': 1 light = 25%, 2 = 50%, 3 = 75%, 4 = 100% battery charge.

All indicator lights off = 100% charged.

## Turning the audible beep on and off

Press and hold both buttons at the same time for 3 seconds. Three short beeps will play;

- High pitch – High pitch – High pitch = audible beep ON.
- High pitch – High pitch – Low pitch = audible beep OFF.

## Warnings

- Do not strap the Device on or near an injury.
- Do not apply any creams or gels to the conductive pads. Only apply Arc conductive gel to the carbon rubber pads, if required.
- Do not apply therapy while driving, operating machinery, or during any activity in which electrical stimulation can put the user at risk of injury.
- The Device is in addition to, and in no way a substitute for, your existing medical treatment or any medical advice you have received from your healthcare professional.
- Always consult a medical expert/healthcare professional on any health conditions and on using the Device device prior to use.
- Like any electrical device, do not use your Device in wet or damp conditions.
- Some people may, after prolonged use, develop a sensitivity reaction. This can manifest in a skin rash type condition. \*\*Therefore, we recommend that the arm being used as well as the positioning of the device is alternated daily.

- Some people are, or can become, sensitised to the conductive pads. In order to minimise the likelihood of any sensitisation, we recommend that the arm being used is alternated daily.
- Device users should stop using the Device and should contact the Arc4Health office if they experience any adverse reactions from the Device. Our contact details are below.
- Please use the Device only as directed.

## Precautions

- Use caution if microcurrent therapy is applied over areas of skin with hypersensitivity.
- Keep Device out of reach of children.

## General Warnings

- Do not immerse any part of the unit in water.
- Use this Device only with the leads, electrodes and accessories recommended by the manufacturer. Use of other parts or materials not supplied by the manufacturer will invalidate the warranty.
- Keep the Device away from fireplaces or radiant heaters as the heat may affect the Device.

- Use of microcurrent therapy may cause some dehydration. Please ensure you drink sufficient water and stay hydrated when using the Device.
- Use of the Device is in addition to, and in no way a substitute for your normal, existing practice of seeking help, guidance and support from your GP, or other Clinical Practitioner.

## Troubleshooting

- Right activation button showing amber = loss of connectivity. Please refer to the Good/Bad connectivity notes below.
- Left On/Off button showing amber = lack of charge. Device unit needs charging. Please refer to the Charging notes above.
- Device keeps switching off = the right activation button was not pressed.
- Audible beep has gone = Device's default setting is to have audible sounds. Please refer to notes for turning the audible beep on and off.

## Good/bad connectivity

If the right play button turns green, the programme has started and you have good connectivity.

If the right activation button turns amber, the programme has started but you DO NOT have connectivity.

If AMBER, please check that;

- The electrical lead from the arm strap is still attached to the bottom of the Device.
- The conductive pads have not come loose or fallen off. The conductive pads have sufficient adhesiveness and have not dried out.

## Specifications

### Arc4Health Delivery Unit

Battery	Li-ion Rechargeable, 3.7V (fully traceable) – non removable
Pad connect	Via 2.1 x 5.5mm (centre positive) plug on unit
Maximum pad voltage	+/- 9V
Programme selection	Via push buttons
Modes	4 programmes Pulse duration: 0.11msec Pulse to pulse spacing: 0.86msec Pulse polarity: bipolar
Output channels	Single channel
Output voltage	0-9V (+/-), controlled to maintain the specified current Max = 4V for R = 10kohm (typical) Max = 0.2V for R = 500 ohms
Output current	0-400 uA
Tolerances	All output parameters less than +/- 10%
Physical dimensions	45mm (W) x 15mm (D) x 105mm (L)

Approximate weight	65g
Operating temperature	0°C to 40°C
Storage temperature	0°C to 55°C
Relative humidity	10 to 65%
IP rating	IP 3X
Expected life of unit	5 years under normal use

### DC power supply

Connection	Via 2.1 x 5.5mm (centre positive) plug on unit (pads removed whilst charging)
Mains supply	100-240 VAC, 50-60Hz, 0.15A
Output supply	12VDC, 0.5A
Compliance	EN 60950

### Warranty

Our limited warranty covers defects in materials and workmanship or electrical fault under normal use for 2 years. This is valid from the date of delivery of the Device to the original purchaser. If a defect arises and a valid claim is received within the warranty period, we will either refund the cost of the defective Device or exchange the Device with a Device that is new or which has been manufactured from new or serviceable used parts.

Exclusions: Conductive Pads, Arm Strap and Charger are not covered by this warranty – only the Device.

This warranty is set out in and subject to our terms and conditions of purchase.

Electromagnetic Compatibility: Precautions

The Device complies with the requirements of EN 60601-1-2, but it is recommended that close proximity to portable and mobile RF communications such devices (WiFi, hubs, mobile and cordless telephones etc) be restricted during treatments as these can influence medical electrical equipment.

### WEEE Statement

The Symbol on the product indicates that it should not be disposed of with your normal household waste but should be taken to a designated collection point for the recycling of waste electrical and electronic equipment.

### Standard Symbols

CE Marking

Manufacturer

Type BF Applied Part

Complies with WEEE regulations

### EC Declaration of Conformity

The Device Microcurrent Therapy System complies with the requirements of the Medical Devices Directive 93/42/EEC as amended. The Device is considered to be a Class IIa Medical Device. Compliance has been demonstrated with reference to the following harmonised standards:

- EN 60601-1
- EN 60601-2-10
- EN 60601-1-2

### Contact details

For further information please direct any general enquiries to:

admin@Arc-Family.com

**Please note:** All enquiries will be handled between 09:00 – 17:00 Monday - Thursday, 09:00 – 16:00 Friday.

Please visit [www.arc4health.com](http://www.arc4health.com) should you need more information on our Device (including our FAQs) or if you need to buy more consumables.

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